

## Web Ambassador™

---

The Web Ambassador is easy to use software to confidentially review patient exams via the Internet. For each exam, you can review DICOM images, a report, and a recorded audio conclusion. The Web Ambassador application can be downloaded to any Internet-connected PC. The Web Ambassador can be used at a physician's office or home, ward, surgical center, clinic, or ER. You will also have the ability to review images, exams, or report summaries using a standard Internet Explorer browser.

Physicians will be able to access all studies (X-Ray, Ultrasound, Mammography, MRI and CT) from Gateway Radiology, Parkway Radiology and Pomerado Imaging on the system.

### Web Ambassador requirements

To install and use Web Ambassador, your personal computer must meet the following requirements:

Item	Requirements
Processor and speed	Pentium-class processor, 1 GHz or higher.
RAM	512 MB (minimum).
Operating system	Microsoft® Windows XP. or Microsoft Windows 2000 with Service Pack 4 or higher .
Internet browser	Internet Explorer 6 or higher.
Available hard disk space	5 GB (minimum) of available hard disk space for program files and to copy exams to the computer's local drive.
Monitor	Video display resolution of at least 1024 x 768 with 256 colors. Note: 15 and 16-bit color modes are not recommended since they do not provide for the display of 256 shades of gray. 24-bit and 32-bit color modes provide for 256 shades of gray and excellent color display, but may have some impact on graphics performance for larger images.
Sound card	Required for audio playback.

**Note:** With many Windows-based applications, performance may increase with additional RAM and processor speed. Contact your system administrator to determine if your computer will gain increased performance from upgrading these items.

